




Product name	Invivo 16 channel Knee Coils
Identification code(s)	4535-301-04181 16 Channel Knee coil 1.5T 4535-301-04191 16 Channel Knee coil 3.0T
Total weight [kg]	6.6 Kg
Manufacturer	Invivo Corporation. USA
Address	3545 SW 47th Avenue
Zip code	Gainesville, FL 32608
Country	USA
Electronic info	RFCoils_Invivo_Knee16ch v00.doc

Recycling information

Special Attention 	Items	Location
	None	Figure

Batteries 	Items, type	Location
	None	Figure

Hazardous substances 	Items	Location
	Several printed circuit boards inside the coil	Figure 1
	Beryllium copper contact pins on connectors	Figure 1,2,3

Fluids and gases 	Items	Location
	Plastic bottle with Nickel Chloride 31.5mM/L of NiCL2-6H2O	Figure 4

Name: Recycling passport Invivo 16channel Knee coils		12NC 4535 301 04180	Version: 00 Doc.Date: 2009-05-21	
	Page 1 of 2	Group No: 5790	Section: -	Status: Authorized

Locations as mentioned in the Passport (Information in pictures)




Figure 1

Figure 2



Figure 3



Figure 4

Name: Recycling passport Invivo 16channel Knee coils		12NC 4535 301 04180	Version: 00 Doc.Date: 2009-05-21
	Page 2 of 2	Group No: 5790	Section: - Status: Authorized