






Printed copies are uncontrolled

<b>Product name:</b>	<b>Bedside CO2/Respiratory Mechanics Monitors</b>	
<b>Identification code(s)</b>	<b>Model 1265, 7100, 8000, 8100</b>	
<b>Total weight (in Kg)</b>	<b>3.63 Kg</b>	
<b>Producer/ Manufacturer</b>	Name company:	Respironics California, INC.
	Address:	2271 Cosmos Court, Carlsbad, CA
	Zip code:	92011
	Country:	United States


<b>Recycle Info</b>	<b>Items:</b>	<b>Location</b>
<b>Special attention</b> 	Disinfect device prior to handling	Figures 1, 2, 3, 4
<b>Fluids / Gases</b> 	<b>Items:</b>	<b>Location</b>
	None	
<b>Batteries</b> 	<b>Type:</b>	<b>Location</b>
	Lead Acid Battery	Figure 6, Location 1
<b>Hazardous</b>  <b>To be Removed</b>	<b>Substances:</b>	<b>Location</b>
	Mercury – Florescent tube backlight contains mercury	Figure 6, Location 5
	PCB – Main board (Lead Solder)	Figure 6, Location 2
	PCB – ETCO2 board (Lead Solder) Models 1265 and 7100 only	Figure 6, Location 3
	PCB – Power supply board (Lead Solder)	Figure 6, Location 4
	PCB – Display board (Lead Solder)	Figure 6, Location 5

Passport Number	Version	Author	ECO	Signature/Date
1016825MC6	A	Patrick Freeman	5-13117	EDMS

Figure 1 Model 1265 CAPNOGARD




Front View




Rear View

Printed copies are uncontrolled

Figure 2 Model 7100 CO<sub>2</sub>SMO




Front View




Rear View

Printed copies are uncontrolled

Figure 3 Model 8000 CO<sub>2</sub>SMO Plus




Front View




Rear View

Printed copies are uncontrolled

Figure 4 Model 8100 CO<sub>2</sub>SMO Plus



Front View



Rear View

Printed copies are uncontrolled

Show locations of materials mentioned on the previous sheet. (Picture information)

Printed copies are uncontrolled




Figure 5 Bottom (All Models)

Printed copies are uncontrolled

