





Recycling Passport

Product name:	Velara GCF 2T	Velara GCF 2T	Velara GFD
Identification code	9890-000-62022	9890-000-62023	9890-000-62091
Product name:	Velara GCF 2T Rel. 3	Velara GCF 2T Rel. 3	Velara GCF 2T R3.9.1
Identification code	9890-000-70301	9890-000-70303	9890-000-70304
Product name:	Velara GCF Rel.2 2T	Velara GCF 2T R.3.11.1	R.3.12.1
Identification code	9890-000-70161	9890-000-70305	
Product name:	Velara GCF 2T R.3.12.1		
Identification code	9890-000-70306, 9890-000-70307		
Order No.			
Device Type	B2B (Business to business)		
Producer /Manufacturer	Company name:	Philips Medical Systems DMC GmbH	
	Address:	Roentgenstrasse 24	
	Zip code:	22335 Hamburg	
	Country:	Germany	
	Web info:	http://www.medical.philips.com/main/company/sustainability/recycling/xray_generators/	
Total weight (in Kg)			315
thereof			

steel, iron	iron, low alloy (<5%)	136,263	KG
	iron, high alloy (>5%)	5,674	KG
steel, iron		141,937	KG
nonferrous metals and alloys	aluminum, -alloy	12,655	KG
	copper, -alloy	5,075	KG
	zinc, -alloy	1,652	KG
nonferrous metals and alloys		19,382	KG
glass / ceramics	glass	0,288	KG
glass / ceramics		0,288	KG
plastics / organic substances	oil	64,06	KG
	thermoplastic	6,524	KG
	thermoset	4,505	KG
	elastomer	0,247	KG
plastics / organic substances		75,336	KG
standard parts	other electronic powered devices	45,92	KG
	printed circuit boards	11,835	KG
	cables	9,173	KG
	power supplies / transformers	6,184	KG
	mounting parts, attaching part	2,17	KG
	fans	1,804	KG
	motors, pumps	0,52	KG
	computers and accessories	0,218	KG
standard parts		77,824	KG
		314,767	KG

Recycling Passport

Recycle info		
General Hazards (Desinfection)	For dismantling activities Treatment Facilities must consider the national requirements. For personnel which can come into contact with contaminated material, preventive measures pursuant to national requirements must be taken into account	
Hazardous Substances  To be removed	Substances:	Location
Batteries  To be removed	Type:	Location
	Lithium chrome cell, 3V	Figure 1 (2)
	Lithium chrome cell, 3V	Figure 1 (3)
Special attention 	Item	Location
	Pb is present in the soldering process of PCBs	Figure 1, 2
Fluids / Gases 	Item	Location
	Transformer Oil, Shell Diala	Figure 1 (1)
	Transformer Oil, Shell Diala	Figure 1 (4)
	HV-transformer oil contains no PCBs	

Recycling Passport




Figure 1