



Designed for high resolution imaging

Philips dStream FootAnkle 8ch coil solution for Ingenia 1.5T / 3.0T

Ski-boot-shaped coil for optimum coverage of the ankle and entire foot down to the toes. The coil design and element layout allow for either large FOV imaging of the whole foot or small FOV high resolution imaging of ankle joints. The coil is easy to set up and can be used with the patient's foot vertical or up to 15 degrees plantar flexed.

- Designed for high resolution visualization of detailed cartilage structures
- Easy patient set up, even for tall patients
- Foot can be positioned in a comfortable 15 degree angle

PHILIPS
sense and simplicity

dS FootAnkle 8ch coil solution

Coverage	30 cm
Max. number of channels	8
Main applications	Foot Ankle Toes
Coil solution type	Integrated
Coil connection	dStream Interface
SENSE parallel imaging	dS-SENSE enhanced parallel imaging performance



Please visit www.philips.com/Ingenia



© 2011 Koninklijke Philips Electronics N.V.
All rights are reserved.

Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

Philips Healthcare is part of Royal Philips Electronics

www.philips.com/healthcare
healthcare@philips.com

Printed in The Netherlands
4522 962 70031 * JAN 2011