



Coverage and high resolution

Philips dStream Knee 16ch coil solution for Ingenia 1.5T

Coil designed for ultra-high SNR imaging over an extended field of view of the knee. Two overlapping rings of eight elements extend the coverage area and minimize the need for precise positioning. dS-SENSE enhanced parallel imaging can be selected in all directions. The dS Knee 16ch coil has a split design for easy patient setup and an ergonomically ramped insert for patient comfort.

- Faster exams with dS-SENSE in all directions
- Large aperture of 15 cm for muscular or obese patients
- Easy positioning and large coverage for high throughput

PHILIPS
sense and simplicity

dS Knee 16ch coil solution

Coverage	20 cm
Max. number of channels	16
Main applications	Knee
Coil solution type	Dedicated
Coil connection	dStream Interface
SENSE parallel imaging	dS-SENSE enhanced parallel imaging performance



Please visit www.philips.com/Ingenia



© 2011 Koninklijke Philips Electronics N.V.
All rights are reserved.

Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

Philips Healthcare is part of Royal Philips Electronics

www.philips.com/healthcare
healthcare@philips.com

Printed in The Netherlands
4522 962 70101 * JAN 2011