



# Designed for high resolution knee imaging

Philips dStream Knee 8ch coil solution for Ingenia 1.5T / 3.0T

Coil designed for high resolution imaging of the left or right knee. The coil closely conforms to the anatomy for optimal signal-to-noise ratio. The coil can be slightly rotated relative to its base plate to ease coil setup and enhance patient comfort.

- Designed for high resolution visualization of fine cartilage structures and small ligaments
- Coil rotation feature for enhanced positioning and comfort
- Split design and single handle for easy coil setup

**PHILIPS**  
sense and simplicity

# dS Knee 8ch coil solution

Coverage	16 cm
Max. number of channels	8
Main applications	Knee
Coil solution type	Dedicated
Coil connection	dStream Interface
SENSE parallel imaging	dS-SENSE enhanced parallel imaging performance



Please visit [www.philips.com/Ingenia](http://www.philips.com/Ingenia)



© 2011 Koninklijke Philips Electronics N.V.  
All rights are reserved.

Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

Philips Healthcare is part of Royal Philips Electronics

[www.philips.com/healthcare](http://www.philips.com/healthcare)  
[healthcare@philips.com](mailto:healthcare@philips.com)

Printed in The Netherlands  
4522 962 70091 \* JAN 2011