



Advancing mainstream body applications

Philips dStream Torso coil for Ingenia 1.5T / 3.0T

An integrated coil solution for body and peripheral vascular related imaging. It includes the FlexCoverage Anterior coil. Combined with the FlexCoverage Posterior coil it allows 60 cm coverage with up to 32 channels.

The flexible, lightweight easy-to-position FlexCoverage Anterior coil is designed to conform both in right-left and foot-head directions for almost any patient. This allows large coverage and comfortable strap-free operation. dS Torso features DirectDigital sampling in the coil, capturing data where the MR signal is at its purest. Digital transmission without losses in the RF chain provides enhanced SNR and dS-SENSE enhanced parallel imaging performance.

- Large 60 cm coverage for ease of patient positioning
- Conforming and strapless design for patient comfort
- Digitization in the coil for increased SNR
- dS-SENSE to speed up acquisition in all imaging directions

PHILIPS
sense and simplicity

dS Torso coil solution

Coverage	60 cm
Max. number of channels	32
Main applications	Torso Chest Pelvis Heart Peripheral-vascular
Coil solution type	Integrated
Data sampling	DirectDigital
Coil connection	Single FlexConnect
SENSE parallel imaging	dS-SENSE enhanced parallel imaging performance in AP, LR and FH directions



Please visit www.philips.com/Ingenia



© 2011 Koninklijke Philips Electronics N.V.
All rights are reserved.

Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

Philips Healthcare is part of Royal Philips Electronics

www.philips.com/healthcare
healthcare@philips.com

Printed in The Netherlands
4522 962 70001 * JAN 2011