



High resolution wrist imaging

Philips dStream Wrist 8ch coil solution for Ingenia 1.5T / 3.0T

Coil that closely fits the left or right wrist for high SNR. This design provides the high SNR needed to acquire images with a small FOV. It has a one piece, ovoid, hinged design for easy patient set up. High quality imaging can be obtained with the coil at the patient's side. The coil attaches to a rigid base plate for fixation to reduce patient motion artifacts.

- Designed for high resolution visualization of interosseous ligaments
- One-piece hinged design for easy patient setup
- Comfortable positioning at the patient's side

PHILIPS
sense and simplicity

dS Wrist 8ch coil solution

Coverage	8 cm
Max. number of channels	8
Main applications	Wrist
Coil solution type	Dedicated
Coil connection	dStream Interface
SENSE parallel imaging	dS-SENSE enhanced parallel imaging performance



Please visit www.philips.com/Ingenia



© 2011 Koninklijke Philips Electronics N.V.
All rights are reserved.

Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

Philips Healthcare is part of Royal Philips Electronics

www.philips.com/healthcare
healthcare@philips.com

Printed in The Netherlands
4522 962 70041 * JAN 2011