



Flexibility for multiple applications

Philips dStream Flex M coil solution for Ingenia 1.5T / 3.0T

An integrated coil solution for general-purpose imaging. It includes two medium-sized flexible general-purpose coils. Combined with the FlexCoverage Posterior coil they allow 15 cm coverage with up to 6 channels.

The shape and size of the flexible coil elements allow a wide variety of applications, including imaging of medium sized anatomies. The coil can be used to locally enhance resolution of images acquired over a larger FOV, for example in pediatric applications.

- Allows positioning close to the anatomy to increase SNR
- Open, non-rigid design offers imaging solutions for challenging cases
- Delivers high resolution and extended coverage in combination with posterior coil

PHILIPS
sense and simplicity

dS Flex M coil solution

Coverage	15 cm
Max. number of channels	6
Main applications	Shoulder Foot Ankle Knee Pediatrics
Coil solution type	Integrated
Coil connection	dStream Interface
SENSE parallel imaging	dS-SENSE enhanced parallel imaging performance



Please visit www.philips.com/Ingenia



© 2011 Koninklijke Philips Electronics N.V.
All rights are reserved.

Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

Philips Healthcare is part of Royal Philips Electronics

www.philips.com/healthcare
healthcare@philips.com

Printed in The Netherlands
4522 962 70071 * JAN 2011