



# Flexibility for high resolution

Philips dStream Small Extremity 8ch coil solution for Ingenia 1.5T / 3.0T

Semi-flexible coil designed for imaging of elbows, hands and small knees. The coil has an inner diameter of 20 cm to match the size of the small extremities. It has a flexible wrap-around design for easy positioning and good fit. A mattress that supports both patient and coil is provided to increase patient comfort and avoid motion.

- Flexible wrap-around design to accommodate various patient sizes
- Close fit to the anatomy enabling high resolution imaging of cartilage and bone
- Positioning pads for patient comfort and immobilization

**PHILIPS**  
sense and simplicity

# dS Small Extremity 8ch coil solution

Coverage	20 cm
Max. number of channels	8
Main applications	Elbows Hands Small knees
Coil solution type	Dedicated
Coil connection	dStream Interface
SENSE parallel imaging	dS-SENSE enhanced parallel imaging performance



Please visit [www.philips.com/Ingenia](http://www.philips.com/Ingenia)



© 2011 Koninklijke Philips Electronics N.V.  
All rights are reserved.

Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

Philips Healthcare is part of Royal Philips Electronics

[www.philips.com/healthcare](http://www.philips.com/healthcare)  
[healthcare@philips.com](mailto:healthcare@philips.com)

Printed in The Netherlands  
4522 962 70051 \* JAN 2011