



Integrated solution for advanced whole body MR

Philips dStream WholeBody coil solution for Ingenia 1.5T / 3.0T

An integrated coil solution for whole body and peripheral vascular related imaging. It includes two FlexCoverage Anterior coils. Combined with the FlexCoverage Posterior, HeadNeck and Base it enables 200 cm coverage with a maximum of 108 channels.

The flexible, lightweight easy-to-position FlexCoverage Anterior coil is designed to conform both in right-left and foot-head directions for almost any patient. This allows large coverage and comfortable strap-free operation.

- Full anatomic coverage up to 200 cm
- Digitization in the coils for increased SNR
- Fast and easy setup with only 3 FlexConnect coil connections

PHILIPS
sense and simplicity

dS Whole-body coil solution

Coverage	200 cm
Max. number of channels	108
Main applications	Whole body Peripheral-vascular Torso Chest Pelvis Heart
Coil solution type	Integrated
Data sampling	DirectDigital
Coil connection	Triple FlexConnect
SENSE parallel imaging	dS-SENSE enhanced parallel imaging performance in AP, LR and FH directions



Please visit www.philips.com/Ingenia



© 2011 Koninklijke Philips Electronics N.V.
All rights are reserved.

Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

Philips Healthcare is part of Royal Philips Electronics

www.philips.com/healthcare
healthcare@philips.com

Printed in The Netherlands
4522 962 70011 * JAN 2011