

## Recycling Passport

<b>Product name:</b>	MammoDiagnost DR		
<b>Identification code</b>	9890-010-87661		
<b>Order No.</b>			
<b>Device Type</b>	B2B (Business to business)		
<b>Producer /Manufacturer</b>	Company name:	Philips Medical Systems DMC GmbH	
	Address:	Roentgenstrasse 24	
	Zip code:	22335 Hamburg	
	Country:	Germany	
	Web info:	See Philips Intranet	
<b>Total weight (in Kg) thereof</b>			<b>508 - 516 Kg</b>

Klasse	Merkmalbezeichnung	Σ	Gesamtgewicht	Eh
steel, iron	iron, low alloy (<5%)		263,580	KG
	iron, high alloy (>5%)		16,710	KG
	cast iron, cast steel		0,447	KG
<b>steel, iron</b>			<b>280,737</b>	<b>KG</b>
nonferrous metals and alloys	aluminum, -alloy		50,922	KG
	copper, -alloy		8,329	KG
	zinc, -alloy		5,120	KG
	nickel, -alloy		0,017	KG
	tungsten, -alloy		0	KG
<b>nonferrous metals and alloys</b>			<b>64,388</b>	<b>KG</b>
precious metals and precious metal alloy	silver, -alloy		0	KG
<b>precious metals and precious metal alloy</b>			<b>0</b>	<b>KG</b>
glass / ceramics	glass		23,266	KG
	aluminum oxide		1,974	KG
<b>glass / ceramics</b>			<b>25,240</b>	<b>KG</b>
plastics / organic substances	thermoplastic		12,799	KG
	thermoset		10,154	KG
	oil		5,117	KG
	elastomer		1,811	KG
	lubricating grease		0,015	KG
<b>plastics / organic substances</b>			<b>29,896</b>	<b>KG</b>
batteries	batteries, lithium		0,008	KG
<b>batteries</b>			<b>0,008</b>	<b>KG</b>
standard parts	bearing plates		49,436	KG
	computers and accessories		17,697	KG
	monitors		16,770	KG
	cables		15,200	KG
	printed circuit boards		6,102	KG
	other electronic powered devic		2,633	KG
	power supplies / transformers		0,601	KG
	mounting parts, attaching part		0,120	KG
<b>standard parts</b>			<b>108,559</b>	<b>KG</b>
relevant materials	lead and -compounds (cables an		5,822	KG
	chromium(VI) -compounds		0,941	KG
	selenium and -compounds		0,070	KG
	cobalt and -compounds		0,005	KG

Printed copies are uncontrolled documents

## Recycling Passport

	antimony and -compounds	0,001	KG
	confirmed by supplier		
relevant materials		6,839	KG
prohibited materials	confirmed by supplier		
prohibited materials			
		515,667	KG


<b>Recycle info</b>		
<b>General Hazards (Desinfection)</b>	For dismantling activities Treatment Facilities must consider the national requirements. For personnel which can come into contact with contaminated material, preventive measures pursuant to national requirements must be taken into account	
<b>Hazardous Substances</b>  <b>To be removed</b>	<b>Substances:</b>	<b>Location</b>
	Lead	Fig. 1 ; 2
	selenium	Fig. 1
	cobolt	Fig. 1
	antimony	Fig. 1
	Chrom (VI)	Fig. 1
<b>Batteries</b>  <b>To be removed</b>	<b>Type:</b>	<b>Location</b>
<b>Special attention</b>  	<b>Item</b>	<b>Location</b>
<b>Fluids / Gases</b>  	<b>Item</b>	<b>Location</b>
	Oil tank / Remove oil before dismantling	Fig. 1
	HV-transformer oil contains no PCBs	

Printed copies are uncontrolled documents

# Recycling Passport

Locations as mentioned in the Passport (Pictures information)


Printed copies are  
uncontrolled documents



MammoDiagnost DR – Basic Unit (508 kg)



- Advanced Stereo option (8 kg)



Console