





PHILIPS HEALTHCARE	Product Recycling Passport	Page 1 of 6
-------------------------------	-----------------------------------	-------------

Product name:	Brilliance CT Gantry	
Identification code(s)	4535 670 23331, 4550 122 04021, 4550 113 09031, 4535 670 79601	
Total weight (in Kg)	2041	
Producer/ Manufacturer	Name company:	Philips Medical Systems
	Address:	595 Miner Road Highland Heights, Ohio
	Zip code:	44143
	Country:	United States of America

Printed copies are uncontrolled

Recycle Info	Items:	Location
Special attention 	Spring loaded cover	Figure 1 (1)
	Actuators: Unstable when removed, possible rotation of frame	Figure 3 (2), Figure 4 (2)
	DMS: System Unbalanced when DMS is removed causing rotation	Figure 2 (6)
Fluids / Gases 	Items:	Location
	Hydrocarbon Oil	Figure 2 (2),(7)
Batteries 	Type:	Location
To be Removed		
Hazardous 	Substances:	Location
	Printed Circuit Boards	Figure 1 (2), Figure 3 (1), Figure 4 (1)
	Lead	Figure 2 (1),(3),(4)
To be Removed		

Title: Recycle Passport Brilliance CT System Gantry
 Recycling Passport Number: 453567149361 Rev: E

Show locations of materials mentioned on the previous sheet. (Picture information)

Printed copies are uncontrolled




Figure 1: Front View of Gantry with Covers on

Printed copies are uncontrolled




Figure 2: Front View of Gantry with Front Cover Lifted

Printed copies are uncontrolled




Figure 3: Right Side View of Gantry with Cover Removed

Printed copies are uncontrolled




Figure 4: Left Side View of Gantry with Cover Removed

Title: Recycle Passport Brilliance CT System Gantry
Recycling Passport Number: 453567149361 Rev: E

Recycling Passport Revision History

Rev	Description	Date
C	CO 2015923 Release of on line documentation for WEEE	2 May 2007
D	CO 2019848 – Changed Identification label to Product Part Number	25 February 2009
E	CO 2020759 – Updated Format	24 June 2009

Printed copies are uncontrolled