





PHILIPS HEALTHCARE	Product Recycling Passport	Page 1 of 3
-------------------------------	-----------------------------------	-------------

Product name:	Brilliance CT Patient Support	
Identification code(s)	4535-670-23322	
Total weight (in Kg)	395	
Producer/ Manufacturer	Name company:	Philips Medical Systems
	Address:	595 Miner Road Highland Heights, Ohio
	Zip code:	44143
	Country:	United States of America

Printed copies are uncontrolled

Recycle Info	Items:	Location
Special attention 	Removal of Motor: Couch will drop when motor is removed!	Figure 1 (1)
	Green Safety Block: When removed, couch may collapse!	Figure 1 (3)
	Giant Spring: Potential energy hazard!	Figure 1 (4)
Fluids / Gases 	Items:	Location
Batteries 	Type:	Location
To be Removed		
Hazardous 	Substances:	Location
To be Removed	Printed Circuit Boards	Figure 1 (2)

Show locations of materials mentioned on the previous sheet. (Picture information)

Printed copies are uncontrolled




Figure 1: Patient Support

PHILIPS HEALTHCARE	Product Recycling Passport	Page 3 of 3
-------------------------------	-----------------------------------	-------------

Recycling Passport Revision History

Rev	Description	Date
C	CO 2015923 Release of Brilliance CT patient support online WEEE passport	2 May 2007
D	CO 2019848 – Updated Identification Code	24 May 2009
E	CO 2020759 - Removed superfluous footer, added titles to supplied pictures, corrected Rev A & B description to be “CO 2019848”	09 June 2009

Printed copies are uncontrolled