

| | | |
|------------|---|-------------|
| Dameca A/S | Dameca MRI 508 Product Recycling Passport | Page 1 of 2 |
|------------|---|-------------|

| | | |
|---------------------------|---|---------------------------|
| Product name: | Dameca MRI508 | |
| Identification code(s) | 10651MRI-00 | |
| Total weight (in Kg) | Approx. 110 kg (depending on configuration) | |
| Producer/ Manufacturer | Name company: | Dameca A/S |
| | Address: | Islevdalvej 211; Roedovre |
| | Zip code: | 2610 |
| | Country: | Denmark |





| Recycle Info | Items: | Location |
|---|---|---|
| Special attention  | Before any maintenance or repair is started the machine must be disconnected from the mains and all gas pressures must be released. Be sure to release the N2O pressure before the O2 pressure. | |
|  Fluids / Gases | Items: Absorber material | Location Figure 1 |
| Batteries  To be Removed | Type: 2 12V Lead acid batteries | Location Figure 1 |
| Hazardous  To be Removed | Substances: PCB boards contain leaded solder: 2 Lead-acid batteries | Location Figure 1 Figure 1 |

Figure 1:

Figure 1:

