





<b>Product name:</b>	<b>Siesta i Whispa</b>		
<b>Identification code(s)</b>	<b>10651-00</b>		
<b>Total weight (in Kg)</b>	Approx. 110 kg (depending on configuration)		
<b>Producer/ Manufacturer</b>	<b>Name company:</b>	<b>Dameca A/S</b>	
	<b>Address:</b>	<b>Islevdalvej 211; Roedovre</b>	
	<b>Zip code:</b>	<b>2610</b>	
	<b>Country:</b>	<b>Denmark</b>	

<b>Recycle Info</b>	<b>Items:</b>	<b>Location</b>
<b>Special attention</b> 	Before any maintenance or repair is started the machine must be disconnected from the mains, and all gas pressures must be released. Be sure to release the N2O pressure before the O2 pressure.	
 <b>Fluids / Gases</b>	<b>Items:</b> Absorber material	<b>Location</b> <b>Figure 1</b>
<b>Batteries</b>  <b>To be Removed</b>	<b>Type:</b> 20-Cell NiMH, 24V, 1.5Ah, (ventilator battery)	<b>Location</b> <b>Figure 2</b>
<b>Hazardous</b>  <b>To be Removed</b>	<b>Substances:</b> PCB boards contain leaded solder	<b>Location</b> <b>Figure 2</b>

**Show locations of materials mentioned on the previous sheet.**

Title: Product Recycling Passport for Siesta i Whispa  
 Recycling Passport Number: 10651 RCP Rev: 50

Copies are uncontrolled

Figure 1:




Figure 2:

