





Product name:	IntelliVue MX400 Patient Monitor, IntelliVue MX450 Patient Monitor IntelliVue MX500 Patient Monitor, IntelliVue MX550 Patient Monitor	
Identification code(s)	866060 (MX400), 866062 (MX450), 866064 (MX500), 866066 (MX550)	
Total weight (in Kg)	4.5kg (MX400), 5.1kg (MX450), 5.2kg (MX500), 6.5kg (MX550)	
Producer/ Manufacturer	Name company:	Philips Medizin Systeme Böblingen GmbH
	Address:	Hewlett-Packard-Str.2
	Zip code:	71034 Böblingen
	Country:	Germany

Recycle Info	Items:	Location
Special attention 		

	Items:	Location
 Fluids / Gases		

	Type:	Location
Batteries  To be Removed	Lithium Ion battery 1x	Figure 2 [1]

	Substances:	Location
Hazardous  To be Removed		

Show locations of materials mentioned on the previous sheet.

(Please use an illustration/drawing of the system, photos are discouraged)

Figure 1





Figure 2

MX400

MX450 / MX500



Recycling Passport Revision History

Rev	Description	Date
A	As issued	30-Jul-2013
B	Added pictures of MX450, MX550	15-Nov-2013
C	Removed template header, reformatted pages	18-Feb-2015