





**Philips Healthcare**

Efficia CMS200 License Keys 453564567571, 453564567581, 453564567591,  
453564567601, 453564409671 Recycling Passport

|                               |   |   |
|-------------------------------|---|---|
| <b>Product name:</b>          | Efficia CMS200 License Key  |   |
| <b>Reference Identifier:</b>  | 453564567571, 453564567581, 453564567591, 453564567601,<br>453564409671 |   |
| <b>Total weight (in Kg):</b>  | < 0.3 kg  |   |
| <b>Producer /Manufacturer</b> | Name company:   | Philips Goldway (Shenzhen) Industrial Inc.            |
|                               | Address:  | No. 2, Tiangong Road,<br>Nanshan District<br>Shenzhen |
|                               | Zip code:   | 518057  |
|                               | Country:  | P.R. China  |

| <b>Recycle info</b>  | <b>Substances:</b>   | <b>Location</b> |
|--|--|-----------------|
| <b>Special attention</b><br> | <b>Indicate any part of the product that could pose a hazard to the recycler and treatment facility during disassembly such as capacitors with residual energy or spring loaded valves and show location</b> |                 |
|  | Battery Door   | Figure 1        |
| <b>Fluids / Gases</b><br>   | <b>Items:</b>  | <b>Location</b> |
|  | <b>Indicate Fluids/Gasses in the product.</b>  |                 |
|  | None   |                 |
| <b>Batteries</b><br>        | <b>Type:</b>   | <b>Location</b> |
|  | <b>Indicate all batteries in the product via battery chemistry and show location</b>   |                 |
|  | Two AA/R6 batteries as back-up power (1.5V Alkaline)   | Figure 1        |
| <b>Hazardous</b><br>        | <b>Item</b>  | <b>Location</b> |
|  | <b>Indicate location of printed circuit boards and mercury in back lights as well as any cadmium, flame retardants (PBB and PBB) and hexavalent chromium in the product and show location</b>                |                 |
|  | None   |                 |

**Additional Information as needed**

| <b>Material content (global)</b> | <b>Weight in Kg</b> | <b>Materials content (continue)</b> | <b>Weight in Kg</b> |
|----------------------------------|---------------------|-------------------------------------|---------------------|
| LG AF342                         | 105 grams           | Two AA/R6                           | 50 grams            |
|                                  |                     |                                     |                     |

**Locations as mentioned in the Passport (Pictures information)**



**Figure 1**  