




Philips Medical Systems / Product Recycling Passport

| | | |
|-------------------------------|--|------------------------------------|
| Product name: | Avalon CTS | |
| Identification code | M2720A (Base Station), M2725A/M2726A/M2727A (Transducer) | |
| Total weight (in Kg) | 2.4 kg (Base Station), 0.11 kg (Transducer) | |
| Producer /Manufacturer | Name company: | Philips Medizin Systeme Boeblingen |
| | Address: | Hewlett-Packard Strasse 2 |
| | Zip code: | Boeblingen 71034 |
| | Country: | Germany |
| | | |

| Recycle info | Substances: | Location |
|---|--------------------|-----------------|
| Hazardous | N/A | |
|  To be removed | | |
| | | |
| | | |
| | | |
| | | |

| Batteries | Type: | Location |
|---|---------------------|-------------------|
|  To be removed | Lithium Ion battery | Figure 3 (item 1) |
| | | |
| | | |
| | | |
| | | |

| Special attention | Item | Location |
|---|-------------|-----------------|
|  | N/A | |
| | | |
| | | |
| | | |
| | | |

| | Item | Location |
|--|-------------|-----------------|
| | | |
| | | |
| | | |
| | | |
| | | |

| Material content (global) | Weight in Kg | |
|----------------------------------|----------------------------------|--|
| | Base Station (M2720A) | Transducer (M2725A, M2726A, M2727A) |
| Battery | N/A | 0.03 |
| Steel | 1.1 | N/A |
| Aluminium | 0.3 | N/A |
| Printed Circuits Boards | 0.6 | 0.04 |
| Plastics etc | 0.4 | 0.04 |
| All other material types | N/A | N/A |

Locations as mentioned in the Passport (Pictures information)



Figure 1




Figure 2




Figure 3

