




Philips Medical Systems / Product Recycling Passport

Product name:	MP40/MP50 Patient Monitor	
Identification code	M8003A/M8004A	
Total weight (in Kg)	6.4 Kg	
Producer /Manufacturer	Name company:	Philips Medizin Systeme Boeblingen
	Address:	Hewlett-Packard Strasse 2
	Zip code:	Boeblingen 71034
	Country:	Germany

Recycle info	Substances:	Location
Hazardous  To be removed	N/A	
Batteries  To be removed	Type:	Location
	Lithium Ion battery (optional: 2x, 1x, 0x)	Figure 2 (item 1)
Special attention 	Item	Location
	Magnesium casting housing.	Figure 2 (item 2)
		Figure 3 (item 2)
		Figure 4 (item 2)
	TFT-Display	Figure 1 (item 3)
	Item	Location

Material content (global)	Weight in Kg	Materials content (continue)	Weight in Kg
Magnesium	1.4	Battery	~1
Aluminium	0.3	Steel	0.3
Plastics etc	0.9	All other material types	~1
Printed Circuits Boards	~1.5		

Locations as mentioned in the Passport (Pictures information)



Figure 1




Figure 2




Figure 4




Figure 3

Screws with TORX 20 head