




Philips Medical Systems / Product Recycling Passport

Product name:	MP 80/MP90 Patient Monitor	
Identification code	M8008A/M8010A	
Total weight (in Kg)	6.4 Kg	
Producer /Manufacturer	Name company:	Philips Medizin Systeme Boeblingen
	Address:	Hewlett-Packard Strasse 2
	Zip code:	Boeblingen 71034
	Country:	Germany

Recycle info	Substances:	Location
Hazardous  To be removed	N/A	
Batteries  To be removed	Type:	Location
	N/A	Figure 1 (item 1)
Special attention 	Item	Location
	Magnesium casting housing-part.	Figure 1 (item 1) Figure 2 (item 1)
	Item	Location

Material content (global)	Weight in Kg	Materials content (continue)	Weight in Kg
Magnesium	1.3	Steel	1.9
Aluminium	1.5		
Plastics etc	0.3		
Printed Circuits Boards	1.4		

Locations as mentioned in the Passport (Pictures information)




Figure 1




Figure 2



Figure 3

1

