




Philips Medical Systems / Product Recycling Passport

| | | |
|-------------------------------|---------------------------|------------------------------------|
| Product name: | M3 Patient Monitor | |
| Identification code | M3046A | |
| Total weight (in Kg) | 3.9 | |
| Producer /Manufacturer | Name company: | Philips Medizin Systeme Boeblingen |
| | Address: | Hewlett-Packard Strasse 2 |
| | Zip code: | Boeblingen 71034 |
| | Country: | Germany |
| | | |

| Recycle info | Substances: | Location |
|---|--------------------|-------------------|
| Hazardous | N/A | |
|  | | |
| | | |
| | | |
| | | |
| | | |
| To be removed | | |
| Batteries | Type: | Location |
|  | NiMH battery | Figure 2 (item 1) |
| | | |
| | | |
| | | |
| | | |
| To be removed | | |
| Special attention | Item | Location |
|  | TFT-Display | Figure 1 (item 2) |
| | | |
| | | |
| | | |
| | | |
| | Item | Location |
| | | |
| | | |
| | | |
| | | |
| | | |

| Material content (global) | Weight in Kg | Materials content (continue) | Weight in Kg |
|----------------------------------|---------------------|-------------------------------------|---------------------|
| Steel | 1.3 | Battery | 0.6 |
| Aluminium | 0.1 | All other material types | 0.2 |
| Plastics etc | 1.1 | | |
| Printed Circuits Boards | 0.6 | | |

Locations as mentioned in the Passport (Pictures information)



Figure 1



2




Figure 2

1

