





PHILIPS HEALTHCARE	Product Recycling Passport	Page 1 of 2
-------------------------------	-----------------------------------	-------------

Product name:	Expression MR200 MRI Patient Monitoring System	
Identification code(s)	Model MR200	
Total weight (in Kg)	33.59 Kg	
Producer/ Manufacturer	Name company:	Philips / Invivo
	Address:	12151 Research Parkway Suite 200 Orlando, Fl
	Zip code:	32826
	Country:	USA

Recycle Info	Items:	Location
Special attention 		

	Items:	Location
 Fluids / Gases		

	Type:	Location
Batteries  To be Removed	989803169491 Batteries (2X) in Base Casting	Figure 1 [1]
	CR Battery in Motherboard	Figure 1 [2]
	9065 Batteries in Wireless Modules (2X)	Figure 2

	Substances:	Location
Hazardous 		

Title: MR200 Recycling Passport
 Recycling Passport Number: A-866120-90041 Rev: A

Copies are uncontrolled

To be Removed		
----------------------	--	--

**Show locations of materials mentioned on the previous sheet.
 (Please use an illustration/drawing of the system, photos are discouraged)**

Figure 1




Figure 2

