





PHILIPS HEALTHCARE	Product Recycling Passport	Page 1 of 4
-----------------------	-----------------------------------	-------------

Product name	Expression Patient Monitor (Model MR400, REF 866185)
Identification code(s)	REF 866185
Total weight (in Kg)	44.75 kg (cart only); 50.2 kg (system)


Producer/ Manufacturer	Name company:	Philips Healthcare
	Address:	3000 Minuteman Rd. Andover, MA
	Zip code:	01810
	Country:	U.S.A

Recycle Info	Items:	Location
Special attention 	Plastics – PC + ABS Blend (total weight = 7.4 kg)	Cart, Lower processing unit, Modules
	Aluminum – Sheet Metal, Die Cast, Sand Cast (total weight = 7.07 kg)	Lower processing unit, Display, Patient Connections
	Zinc – Die Cast (total weight = 15.79 kg)	Base Casting
	Glass – Fiber Optic, Touch Screen (total weight = 0.82 kg)	Display, Cart
	Copper – Cables (total weight = 0.28 kg)	Lower processing unit, Display, Patient Connections
	Printed Circuit Boards (total weight = 1.12 kg)	Lower processing unit, Display, Patient Connections, Modules
	Packaging: - PE Foam – 3.65 kg - Wood – 12.98 kg - Cardboard – 6.04 kg	Packaging


 Fluids / Gases	Items:	Location
	N/A	

 <p>Batteries</p> <p>To be Removed</p>	Type:	Location
	989803169491 - (Panasonic) CGA103450A Lithium Ion (12 cell per pack; 4 packs per system) Weight per pack = 637.5 g	Cart lower processing unit
	989803191341 - (BAK) 383562PN Lithium Polymer Cell/Battery (1 cell per pack; 3 packs per system) Weight per pack = 30.2 g	Modules (2 packs for ECG, 1 pack for SpO2)
<p>Hazardous</p>  <p>To be Removed</p>	Substances:	Location
	PVC	Lower processing unit internal tubing
	PUR: - 989803183221 - 989803183231	Patient NIBP Trunk Hoses

Show locations of materials mentioned on the previous sheet.



Copies are uncontrolled



Recycling Passport Revision History

Rev	Description	Date
A	As Issued	2015 May 07

Copies are uncontrolled