





Philips Medical Systems/Product Recycling Passport

Printed copies are uncontrolled

Product name:	PWD Transmitter	
Identification code(s)	862108	
Total weight (in Kg)	0.2 Kg	
Producer/ Manufacturer	Name company:	Philips Medical Systems
	Address:	3000 Minuteman Road, Andover, MA
	Zip code:	01810
	Country:	USA

Recycle Info	Items:	Location
Special attention 	Not applicable	
Fluids / Gases 	Not Applicable	
Batteries 	Type:	Location
	2 AA Alkaline Batteries (optional)	See figure 4
To be Removed		
Hazardous 	Substances:	Location
	Radio PC Board	See figure 3
	Power Supply PC Board	See figure 3
	System PC Board	See figure 2
	Sensor PC Board (optional)	See figure 1
To be Removed		

Title: Transmitter Recycling Passport
 Recycling Passport Number: 453564044681 Rev: A

Philips Medical Systems/Product Recycling Passport

Printed copies are uncontrolled

Figure 1. Sensor PC Board (optional)




Figure 2. System PC Board





Figure 3. Radio and Power Supply Boards



Radio Board

Figure 4. Batteries (optional)



Title: Transmitter Recycling Passport

Recycling Passport Number: 453564044681 Rev: A