





Philips Medical Systems/Product Recycling Passport

Printed copies are uncontrolled

Product name:	Philips Telemetry System Frequency Converter	
Identification code(s)	M2616A	
Total weight (in Kg)	1.8 Kg	
Producer/ Manufacturer	Name company:	Philips Medical Systems
	Address:	3000 Minuteman Road, Andover,MA
	Zip code:	01810
	Country:	USA

Recycle Info	Items:	Location
Special attention 	Not applicable	
Fluids / Gases 	Not Applicable	
Batteries 	Not Applicable	
To be Removed		
Hazardous 	Substances:	Location
	Frequency Converter PC Board	See Figure 1
	Power Supply PC Board	
To be Removed		

Title: Tuhfman External Frequency Converter Recycling Passport
 Recycling Passport Number: M2616-90016 Rev: A

Philips Medical Systems/Product Recycling Passport

Figure 1 PC Boards



Power Supply PC Board

Frequency Converter PC Module

Printed copies are uncontrolled