

---

**Philips Healthcare**

Product Speaker Bar 453564381621 for Philips SureSigns Central 863291  
Recycling Passport

---

**Notice**                      ***UNCONTROLLED IN PRINTED FORM UNLESS STAMPED IN RED***





---

**Revision Summary**            The current document revision is described below:

Rev.	Change Summary    (The CO number of this revision is shown in the DMS web 'effective area' screen for this document)
A	As released

**Document Components**    This document contains no additional components in DMS.

<b>Product name:</b>	<i>Enter name of product as it appears on the product.</i> 453564381621 Speaker Bar	
<b>Reference Identifier:</b>	<i>This should be the model number as it appears on the product.</i> 453564381621	
<b>Total weight (in Kg):</b>	0.354	
<b>Producer /Manufacturer</b>	Name company:	Hewlett Packard
	Address:	3000 Hanover Street, Palo Alto, CA
	Zip code:	94304-1185
	Country:	USA

<b>Recycle info</b>	<b>Substances:</b>	<b>Location</b>
<b>Special attention</b> 	<b>Indicate any part of the product that could pose a hazard to the recycler and treatment facility during disassembly such as capacitors with residual energy or spring loaded valves and show location</b> None	
<b>Fluids / Gases</b> 	<b>Items:</b>	<b>Location</b>
	None	
<b>Batteries</b> 	<b>Type:</b>	<b>Location</b>
	<b>Indicate all batteries in the product via battery chemistry and show location</b> None	
<b>Hazardous</b> 	<b>Item</b>	<b>Location</b>
	<b>Indicate location of printed circuit boards and mercury in back lights as well as any cadmium, flame retardants (PBB and PBB) and hexavalent chromium in the product and show location</b> There are two PCBs	Figure 1

#### Additional Information as needed

<b>Material content (global)</b>	<b>Weight in Kg</b>	<b>Materials content (continue)</b>	<b>Weight in Kg</b>
Lead	<0.1 %		
Plastic	0.181		
Metal	0.136		
Printed Circuits Boards	0.032		

Locations as mentioned in the Passport (Pictures information)



Figure 1

