

POWER TO SAVE A LIFE

HEARTSTART
DEFIBRILLATORS



HEARTSTART TRAINER

TRAINING FOR LIFE

PHILIPS



**PHILIPS' HEARTSTART TRAINER PROVIDES
A REALISTIC TRAINING EXPERIENCE FOR
USERS OF THE PHILIPS HEARTSTART ONSITE
AND HEARTSTART HOME DEFIBRILLATORS.**

Empower your students by teaching them how to use HeartStart Defibrillators to save lives. The HeartStart Trainer simulates how the HeartStart OnSite and HeartStart Home Defibrillators will operate during a real-life sudden cardiac arrest event. Clear, confident

voice instructions that provide easy-to-follow directions, and preconfigured scenarios, make training students easier and more rewarding. And for safety, the HeartStart Trainer is sufficiently different to avoid inadvertent use in an emergency.

SPECIFICATIONS

Standard Features

- ▶ Voice instructions match those of the HeartStart OnSite and HeartStart Home Defibrillators.
- ▶ Works with internal and external pads adapters, allowing the instructor to provide some or no clues toward pad placement; resulting in the most realistic training scenario.
- ▶ Resembles the HeartStart OnSite and HeartStart Home Defibrillators for a realistic training experience, yet is sufficiently different to avoid inadvertent use in an actual emergency.
- ▶ Accepts Infant/Child Training Pads for pediatric defibrillation training.
- ▶ CPR coaching helps train students in adult and infant/child CPR.
- ▶ Battery well and real LED status indicator for maintenance training.
- ▶ Uses four standard AA cell batteries (not included) allowing you to manage training costs.

HEARTSTART TRAINER

M5085A

Comes with nylon carrying case, one set of Adult Training Pads and one External Manikin Adapter Strip for use with standard manikins.

Dimensions:

2.8" H x 8.3" W x 8.3" D

Weight:

2 lbs with batteries

Batteries:

4 AA cells

Training scripts:

Eight preconfigured scripts based on common training scenarios.

Replacement Items and Accessories

REPLACEMENT CARRYING CASE

M5087A

EXTERNAL MANIKIN ADAPTER FOR USE WITH ALL MANIKINS (10 PACK)

M5089A

INTERNAL MANIKIN ADAPTER FOR USE WITH LITTLE ANNE™ AND RESUSCI® ANNE CPR MANIKINS

M5088A



Training Pads

- ▶ Reusable.
- ▶ Resemble defibrillation pads for a realistic training experience, yet sufficiently different to avoid inadvertent use in an emergency.
- ▶ Easy switching between adult and infant/child pads cartridges.

ADULT TRAINING PADS CARTRIDGE FOR HEARTSTART TRAINER

M5073A

INFANT/CHILD TRAINING PADS CARTRIDGE FOR HEARTSTART TRAINER

M5074A



REPLACEMENT ADULT TRAINING PADS, WITH REPLACEABLE SEAL

M5093A

REPLACEMENT INFANT/CHILD TRAINING PADS, WITH REPLACEABLE SEAL

M5094A

Philips, Leading the Way

Philips, the maker of HeartStart Defibrillators, is the global leader and trusted choice for portable defibrillation technology solutions.* Our complete family of resuscitation products helps responders treat sudden cardiac arrest quickly and effectively wherever it occurs—anywhere people live, work and play.

Why Philips is the trusted choice:

- ▶ One of the world's largest medical products companies with annual revenue of over \$6 billion.
- ▶ World market-share leader in automated external defibrillators (AEDs)*
- ▶ First to develop an automated external defibrillator solution to defibrillate infants and children.
- ▶ More than 120,000 Philips defibrillators deployed on airlines, and in airports, workplaces, buildings, schools and communities worldwide.
- ▶ Winner of *BusinessWeek/USA Today's* "Best Product of 2002."
- ▶ Over 2 billion HeartStart Defibrillator service hours logged, with an additional 2.7 million hours added every day.
- ▶ The SMART Biphasic Waveform is backed by more published manuscripts demonstrating safety and efficacy than all other brands combined.
- ▶ Advanced the out-of-hospital early defibrillation movement with the launch of its first automated external defibrillator in 1996, establishing the standard for easy-to-use AEDs.

Learn how the Philips HeartStart Trainer can enhance the professional training you provide by calling 1-800-453-6860 or visit www.medical.philips.com/goto/heartstart

* Frost and Sullivan

PHILIPS MEDICAL SYSTEMS IS PART OF ROYAL PHILIPS ELECTRONICS

On the Web: www.medical.philips.com

Via e-mail: medical@philips.com

By postal service:

Philips Medical Systems
 Global Information Center
 I.B.R.S. / C.C.R.I. Numéro 11088
 5600 VC Eindhoven
 Pays-Bas / The Netherlands

CORPORATE HEADQUARTERS

Philips Medical Systems
 Nederland B.V.
 Postbus 10.000
 5680 DA Best
 Netherlands

UNITED STATES

Andover, Massachusetts
 (800) 934-7372

CANADA

Markham, Ontario
 (800) 291-6743

EUROPE, MIDDLE EAST AND AFRICA

Böblingen, Germany
 Fax: (49) 7031 463 1552

LATIN AMERICA HEADQUARTERS

Sunrise, Florida
 (954) 835-2600

ASIA PACIFIC HEADQUARTERS

Wanchai, Hong Kong
 (852) 2821 5888

©Philips Medical Systems

North America Corporation, 2003

All rights reserved. Reproduction in whole or in part is prohibited without the prior written consent of the copyright holder.

Philips Medical Systems North America Corporation reserves the right to make changes in specifications or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

Resusci Anne and Little Anne are either trademarks or registered trademarks of Laerdal Medical Corporation, Stavanger, Norway.

Printed in U.S. (May 2003)
 5990-0779EN