PHILIPS

Ultrasound

Compact 5000 Series



MSK and procedures

Compact without compromise

Philips Compact Ultrasound System 5000 Series for MSK and procedure POC

Some compact ultrasound systems have capabilities as small as their size. Philips Compact 5000 Series for point of care is different. It's compact without compromise, supporting the procedural guidance and assessment needs of busy musculoskeletal, vascular, breast and endocrine or interventional ultrasound practices.

Advanced performance for assessment and procedures

Compact 5000 Series includes many of the advances of Philips premium ultrasound–including a vast array of Philips high-frequency linear transducers such as the PureWave eL18-4, L12-3ERGO and L15-7io – in a portable system with proven workflow for valuable insights at the point of care.

Portability

Bring advanced scanning capabilities such as MicroFlow Imaging (MFI) and the eL18-4 transducer with needle visualization for assessment and procedural guidance with a compact system that has shared architecture with Philips premium ultrasound systems.

Performance

Enhance your ability to reach a definitive diagnosis or have clear guidance during procedures with Philips advanced ultrasound capabilities, imaging algorithms and cutting-edge technologies such as PureWave transducers, xRes and SonoCT for exceptional image quality even in technically difficult patients.

Proven workflow

Capitalize on the shared architecture of Philips EPIQ and Affiniti ultrasound systems to potentially scan more patients, more quickly with familiar user interfaces, touchscreens, features and workflow, offering high productivity with premium-level clarity.

Bring care where needed

Experience exceptional image quality through the PureWave eL18-4 transducer with needle visualization.

Scan comfortably with favorable system ergonomics, including the 10" touchscreen with adjustable viewing angle.

Tap into wider expertise directly from the ultrasound system using Collaboration Live with multi-party* tele-ultrasound for access to other team members for consultation and training.

Simplify service and maintenance with a modular system with familiar architecture and transducers.

Compact 5000 Series features award-winning mobility.






Quickly boot up from transport mode and scan for up to 2 hours with extended batteries.

Support your infection-control practices with a sealed control panel that is responsive to a gloved touch and is compatible with a range of disinfectants and cleaning solutions.†

Easily get your data where it needs to go with DICOM structured reporting.

Choose a flexible configuration of laptop with cart and/or travel case.‡

Compatible with transducers across the Philips portfolio to make the most of your investment



* Contract required. Diagnostic use, remote access via mobile device or browser, multi-party feature and system-to-system connect require release 2.0 or higher. Not intended to be used with QApps.
† https://www.philips.com.au/healthcare/resources/feature-detail/ultrasound-care-and-cleaning.
‡ Travel cases sold separately.

© 2025 Koninklijke Philips N.V. All rights are reserved.
Philips reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

www.philips.com
Printed in the Netherlands.

00000908-00-01 * MAR 2025